Caddis Puppa

Hook	TMC2487 10 12
	Kamasan B100 10 12
	JMC CA40 10 12
Thread	Schwarz 06 08
Body	Leather strip 3-4mm
	2 Ostrich Herls white
Wing	CDC light brown or olive
	Hare mask brown






1. Wind 7-8 turns of lead around the hook shank



3. Tie in the long Ostrich herls on each side of the hook shank



2. Wind the thread from the eye of the hook into the curvature of the hook



4. Tie in the leather strip at the lead wire an fix it down into the curvature of the hook



- 5. wind the leather strip for half a turn until it points to you
- 6. Wind the Ostrich herl which is on the hook side facing you - around the leather strip (counterclockwise) until its tip points back to the hook curvature
- 7. wind the leather strip for half a turn until it points away from you (this secures the Ostrich herl on your side)
- 8. On the hook side which is opposite to you, we wrap the second Ostrich herl around the leather strip (clockwise) until its tip points back to the hook curvature
- 9. Repeat steps 5 to 8 up to 5mm in front of the eye of the hook.
- 10. Secure the leather strip and cut off the excess (leather and Ostrich herls)



- 12. Insert the CDC fibers into a dubbing loop or split thread
- 13. twist
- 14. wrap around head to form a kind of hackle



11. Use a black permanent marker to color the back of the nymph.



- 15. Insert the hare's ear hair into a dubbing loop or split thread
- 16. twist
- 17. wrap around head to form a kind of hackle
- 18. wind head
- 19. whipfinnish
- 20. Vernis